

Mechanical ventilator

Country of origin | United States of America

Primary function | Treatment

Category | Medical device

Commercial information

List price (USD): 4000

Year of commercialization: N/A

Number of units distributed: N/A

Currently marketed in: N/A

Brand: OneBreath, Inc.



Product description

OneBreath's ventilator addresses the need for continuous ventilation in the absence of compressed gas or electricity and provides a range of settings utilizing a portable O₂ concentrator. Preventative maintenance is minimal and a simple and intuitive user interface controls state-of-the-art clinical features including Volume Targeted Pressure Control (VTPC), spontaneous breath support modes, programmable PEEP, neonatal ventilation, and accurate blending and sensing capabilities.

Product details

Accessories: Includes AC power cord and oxygen sensor

Consumables: The disposable breathing circuits, endotracheal tubes, bacterial filters (optional)

Lifetime: 2-5 years

Energy requirements: Rechargeable battery, Continuous power supply, AC, 110V, 220V, 4-hour battery recharge cycle, 24-hour battery life

Contact: Matthew Callaghan | **Email:** info@onebreathventilators.com | **Phone:** +1 650 646 4984 |

Web: www.onebreathventilators.com

NOTE: Information reported by manufacturer before 17 December 2021

WHO ASSESSMENT

Clinical assessment

In many resource-limited settings, lower respiratory tract infections, along with COVID-19 and respiratory diseases, are a major source of morbidity and mortality, ranking among the top 10 causes of death in low- and lower-middle-income countries. The COVID-19 pandemic has heightened the global demand for high-quality, safe, and affordable medical equipment to treat patients with respiratory failure, particularly those who require mechanical ventilation. Pending regulatory and technical evaluation, OneBreath may provide a suitable solution for mechanical ventilation for patients with respiratory failure, provided that the healthcare workforce receives appropriate, relevant training in its clinical use and maintenance.

WHO specification comparison







OneBreath ventilator partially complies with "Transport ventilator", "Sub-acute ventilator" and "ICU ventilator" WHO technical specifications.

The following main non-compliances and non-verifiable (or not specified) requirements are listed:

























Non-compliance: Volume Control Ventilation (VCV) mode is required. Real-time scalar waveform for volume not available. Tidal Volume expired not displayed. ONLY FOR "ICU VENTILATORS": Tidal Volume range required at least 20-1500 ml, while only 20-1000 ml is available. Three scalar waveforms are required to be displayed at the same time while only two are simultaneously available. Loop (axis) displays for pressure-volume and flow-volume not available. Expiration time and End-Tidal CO₂ parameters not displayed

Not possible to be verified: Oxygen-air mixture accuracy. Internal function testing/leak testing. "IP" protection level certification. Display brightness and contrast control capabilities. Not clear if "Inspiratory pause manoeuvre capability to measure plateau pressure" is available. Not clear if "PEAK, PLATEAU and MEAN airway pressure" and "RR (spontaneous and mechanical)" parameters can be displayed and if "adjustable peak pressure limitation/pressure-cycling mechanism above measured peak pressure" is available. "High/Low FiO₂", "High/Low RR", "High/Low Tidal Volume" and "Low PEEP" alarms availability. Availability of single-limb breathing circuits for adult/paediatric.

Regulatory assessment

































 Pre-market assessment	 Not acceptable	Adequate documentation was not provided to perform a medical device Regulatory or Quality System review on this product. Design verification and validation reports were not submitted for pre-market assessment. No documentation available on quality management system, ISO 13485:2016 for quality system assessment. No documentation was provided on post-market activities. As such, this product is in the prototype stage.
 Post-market assessment	 Not acceptable	
 Quality system assessment	 Not acceptable	

Technology evidence assessment

Domains	Evidence assessment		
	Risk/benefit ratio	Impact	
 Medical			The device is a patent-pending prototype made in the United States and Singapore. There is no provision for local manufacture. According to manufacturer tests, the device is simple to use and maintain in rural locations. It is portable and can be used in multiple settings, including ICUs, Emergency rooms, and transportation. Despite this, the device has yet to receive regulatory approval, and the manufacturer has not provided any clinical study data. The unit costs appear to be reasonable, yet they are not inexpensive. In this aspect, affordability is unlikely to be guaranteed in LMI settings. The device is a prototype. In conclusion, a recommendation could not be made based on the current level of evidence.
 Safety			
 Economy			
 Organizational			
 Legal			
 Social			
 Ethical			
 Green environment			

Summary	
Innovation	Technology evidence assessment
Technology readiness level	Not recommended, still a prototype
5	

Health technology and engineering management

Domains	Appropriateness	Domains	Appropriateness	Target settings:
 Durability		 Ease of cleaning		Rural, Urban, Indoors, Ambulance, Secondary & Tertiary level
 Ease of Use		 Ease of maintenance		
 Positive impact on clinical outcomes		 Infrastructure requirements		The core design is based on proven innovation (turbine). The user interface is intuitive. Internal disinfection of the ventilator between patients' use is not addressed. The manufacturer mentions periodic maintenance, however, there is no information on training or parts availability in the submission. Disposables that are not proprietary are indicated. A field evaluation is pending.
 Affordability		 Local access to sales support		
 Engineering resources minimization		 Local access to technical support		
 Cultural and social acceptability		 Local access to training		
 Environmental conditions		 Local access to spare parts		
 Aesthetics		 Locations of use within target setting		

Intellectual property and local production

 **Technology transferability**

 **Openly access intellectual property**

 **Local production**



Intellectual property - It is patent-protected. Some patents are owned by third parties. All applicable patents will require clearance before they can be used. The breath delivery and ventilator management code is proprietary. For the user interface, third-party authorisation is also necessary. Proceed with due diligence.

Local production - It is not yet ready for production; it is still in the prototype stage.

WHO related guidance material

- WHO Global Health Estimates (the top 10 causes of death) - <https://www.who.int/news-room/fact-sheets/detail/the-top-10-causes-of-death>
- WHO Coronavirus (COVID-19) Dashboard - <https://covid19.who.int/>
- Coronavirus disease (COVID-19) technical guidance: Patient management - <https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/patient-management>
- Living guidance for clinical management of COVID-19 - <https://www.who.int/publications/i/item/WHO-2019-nCoV-clinical-2021-2>
- Therapeutics and COVID-19: living guideline - <https://apps.who.int/iris/handle/10665/345356>
- WHO-ICRC Basic Emergency Care: approach to the acutely ill and injured - <https://www.who.int/publications/i/item/basic-emergency-care-approach-to-the-acutely-ill-and-injured>
- Emergency Care - <https://www.who.int/emergencycare/systems/en/>
- Guidelines for essential trauma care - <https://www.who.int/publications/i/item/guidelines-for-essential-trauma-care>
- Severe Acute Respiratory Infections Treatment Centre - <https://www.who.int/publications/i/item/10665-331603>